(English version on page 2)

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Samples

- ⇒ No sample representing a risk for the manipulators or for the equipment will be accepted (radio-marked sample, blood, etc.). In case of doubt, the User must inform the PGTB.
- ⇒ Dilute DNA at ~5-15 ng/μL after quantification by absorbance or fluorescence. If the concentration is < 5 ng/μL or > 15 ng/μL, please notify your PGTB contact by email. Caution: above 30 ng/μL, the quality of the results may be affected.
- ⇒ Dispense 30 μL of the diluted DNA into the following 96 plate references, leaving H6 and H12 wells on each plate empty (without H₂O or DNA) for the PGTB controls, and seal the plates according to the mode of transport that will be used.

Format	Producer	References
96 Plates	4TITUDE	4TI-0770/C, 4TI-0900/C
	THERMOFISCHER	AB-0800, AB-1400, AB-2400
Films and caps	4TITUDE	4Ti-0500, 4TI-0751
	THERMOFISCHER	AB-0558
	EPPENDORF	0030124847

Mode of transport	Sealing of 96 plates
By plane	Cap strips + Parafilm
Road transport	Strongly adhesive film

- ⇒ Identify the plates using the project acronym and the plate number (PROJECT - PLATE 1).

Sample plate information

- ⇒ Use only unique ID for samples, controls and replicates, and use only letters, numbers and underscores, without space to identify samples.
- ⇒ Prepare the sample plate information following this example:

	1	2	3	4	5	6	7	8	9	10	11	12
A	E_32	JK										
B	P_46	LMTR										
C	L_12	AFR										
D	XX											
E	P_89											
F	09											
G	67											
H	AB					empty						empty

- ⇒ Send this file to your contact at PGTB by email.

NEVER SEND STOCK DNA, ONLY ALIQUOTS