(English version on page 2)

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
### Sample preparation

- ⇒ No sample representing a risk for the manipulators or for the equipment will be accepted (radio-marked sample, blood, etc.). In case of doubt, the User must inform the PGTB.
- ⇒ After making your PCR, check on an electrophoresis gel that it is correctly amplified to the expected sizes. Annotate your gel as shown here by specifying:
  - type of size marker used
  - the name of the samples or their position on the plate
  - the expected fragment size
- ⇒ Dispense (at minimal) 25 µL of the product PCR into the following 96 plate references, leaving the H12 well on each plate empty (without H<sub>2</sub>O or DNA) for the PGTB controls, and seal the plates according to the mode of transport that will be used.

Format	Producer	References
96 Plates	4TITUDE	4TI-0770/C, 4TI-0900/C
	THERMOFISCHER	AB-0800, AB-1400, AB-2400
Films and caps	4TITUDE	4Ti-0500, 4TI-0751
	THERMOFISCHER	AB-0558
	EPPENDORF	0030124847

Mode of transport	Sealing of 96 plates
By plane	Cap strips + Parafilm
Road transport	Strongly adhesive film / heat-sealed film

- ⇒ Identify the plates using the project acronym and the plate number (PROJECT - PLATE 1).

### Sample plate information

- ⇒ Use only unique ID for samples, controls and replicates, and use only letters, numbers and underscores, without space to identify samples.
- ⇒ Prepare the sample plate information following this example:

	1	2	3	4	5	6	7	8	9	10	11	12
A	E_32	JK										
B	P_46	LMTR										
C	L_12	AFR										
D	XX											
E	P_89											
F	09											
G	67											
H	AB											empty

- ⇒ Send this file and your annotated gel to your contact at PGTB by email.