(English version on page 2)

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**Pool preparation**

- ⇒ No sample representing a risk for the manipulators or for the equipment will be accepted (radio-marked sample, blood, etc.). In case of doubt, the user must inform the PGTB.
- ⇒ After making your pool, check **its concentration** by qPCR ideally, if not by fluorescence, as well as **its size**.
- ⇒ Place an aliquot (minimum) of **30 µL** of the **5nM** pool in a 1.5mL or 2mL low-bind Eppendorf tube, close and seal the tube with parafilm.
- ⇒ Identify the tube with a label bearing the project acronym and a unique number if there are several tubes.

**NEVER SEND STOCK POOL, ONLY ALIQUOTS**

**Sample information in the pool**

- ⇒ Use only unique identifiers for your samples, controls and repeats, and use only letters, numbers and underscores, with no spaces or dots. Identifiers must not exceed 50 characters.
- ⇒ Prepare a file containing the unique sample ID and indexes (1 index = 1 sample) corresponding to each sample, specifying the desired number of cycles as well as the library loading concentration and the theoretical pHIX percentage :

- *Read1Cycles* : 50
- *Index1Cycles* : 8
- *Index2Cycles* : 8
- *Read1Cycles* : 50

*Loading concentration* : 750pM  
*pHiX percentage* : 5%

In a **single** sequencing run, indexes must ALL be of the same size and format (dual/single).

| Sample ID | UDI            | Index i7 | Index i5 |
|-----------|----------------|----------|----------|
| E_32      | <b>UDI0001</b> | AATCGTTA | ACGTTATT |
| P_46      | <b>UDI0013</b> | AATCGGCG | TCAATATT |
| L_12      | <b>UDI0025</b> | AATCCGTT | AGCATATT |
| XX        | <b>UDI0037</b> | AAGATACA | CAGATATT |
| P_89      | <b>UDI0049</b> | AAGACGAA | AATCTATT |
| 09        | <b>UDI0061</b> | AAGTAAGT | GGACTATT |
| 67        | <b>UDI0073</b> | AAGTTATC | TTGCTATT |

**Note:** If you do not specify a loading concentration, we will load the cassette according to Illumina recommendations with a library concentration of 750pM and a pHIX concentration of 5%.

- ⇒ Send this file with sample <-> index matches to [pgtb-illu@inrae.fr](mailto:pgtb-illu@inrae.fr) by e-mail, together with the **concentration** and **average pool size**, and the **references of the library and index kits** used to prepare the libraries.